



BRN Radial Band



BRN Radial Band



- Unique design prevent unintended movement of the bandage.
- Soft material fits well to patient's skin while easy to clean for the outer layer.
- Double balloon design to achieve better hemostatic effect.
- Puncture site visibility provided by transparent material.
- Standard luer lock connector, compatible with normal syringe.
- Two length available: 24cm, 27cm.
- Include 20ML Syringe.

Ref. NO.	Product Description
V-BRTB24	BRN Transradial Band 240mm
V-BRTB27	BRN Transradial Band 270mm



HOSMED-INC

3403 NW 82nd Avenue. Suite #102.

Doral, Florida 33122. USA

www.hosmed-inc.com

P: (954) 505-1432

E-mail: marketing@hosmed-inc.com

H O S M E D