

“Results matter.”



VITROS® 350

Chemistry System



Designed to help you simply and effectively deliver the results that matter to the patients who need them.

A Compact, Powerful System that Delivers Big Results

The VITROS® 350 Chemistry System offers a range of important benefits such as:

Advanced Security

- + Unique sample processing capabilities, including single use disposable tips, as well as bubble and clot detection for worry-free result reporting
- + MicroSlide technology minimizes interferences by introducing masking and scavengers layers to capture interfering substances

Optimized Efficiency

- + More than 60 onboard reagent positions with a wide accessible menu to deliver immediate results for STAT or routine testing
- + Walk-away operation and simple procedures that allow around-the-clock availability to all members of the lab

Cost Effectiveness

- + Low cost per reportable result and high reagent efficiency without the costs, maintenance, preparation, carryover and interference associated with traditional water-based systems
- + Standardized reagents use the same reagent cartridges across all VITROS® Systems

Key Features

- + Proven VITROS® MicroSlide technology reduces interference
- + Short sample, clot, and bubble detection offers greater accuracy
- + Single use disposable tips eliminate carryover and cross contamination
- + Continuous process verification helps ensure correct results the first time
- + MicroSlide technology and onboard dilution allow random testing of serum, plasma, and urine at any time
- + No plumbing, drains, vents or deionized water eliminates many costs associated with alternative systems
- + Low total costs with increased reportable result efficiency, fewer errors and interventions, and more effective staff utilization
- + Small footprint requires less space and overhead expenses
- + “Load-and-go” reagent preparation, MicroSlide technology and calibration stability up to lot change mean labor can be dedicated to value-added tasks
- + Easy to train and cross-train operators
- + Throughput up to 300 results per hour*

* Depending on test mix



Intellicheck®



MicroSlide

Acetaminophen	Phosphorus
Albumin	Potassium
Alcohol	Salicylate
ALK PHOS	Sodium
ALT	Theophylline
Ammonia	TIBC
Amylase	Total Bilirubin
AST	Total Protein
BuBc Bilirubin	Triglycerides
BUN	Uric Acid
Calcium	Urine Protein
Carbamazepine	
CO2	
Chloride	
Cholesterol	
Cholinesterase	
CK	
CK-MB	
Creatinine	
CRP	
CSF Protein	
dHDL	
Digoxin	
GGT	
Glucose	
Iron	
Lactate	
LDH	
Lipase	
Lithium	
Magnesium	
Phenobarbital	
Phenytoin	

Product availability subject to local regulatory requirements.

Ortho Clinical Diagnostics

PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

VITROS® ^{System} 350
Chemistry

© Ortho-Clinical Diagnostics, Inc. 2013 CL12202
All trademarks are the property of Ortho-Clinical Diagnostics, Inc.